

MICROTITER STAINING METHOD

Use "V" rigid plates, microtiter carriers, IEC Head # 976 with UK centrifuge or 276 head with PR-2 or PR-J IEC.

1×10^6 or 3×10^6 cells may be stained per well (For 1×10^6 add 200λ per well at $5 \times 10^6/\text{ml}$; for 3×10^6 add 200λ per well of cells at $1.5 \times 10^7/\text{ml}$.)

Balance microtiter plates; centrifuge at 1200 rpm ($\approx 290g$) for 3 min. (Using PR-2 upstairs, set rheostat on 27 when using 2 carriers)

timer for 6 minutes
Flick with wrist action to remove supernatant; invert plate and allow to drain on paper towel

Resuspend pellet before adding any liquid. If using microtiter mixer, set on #4 for 3 min. Hand held vibrator may also be used by attaching vibrator to one hand holding plate ~~with the other~~ with the other. Cells may also be resuspended in 50λ liquid by blowing at the pellet in the center of microtiter well above surface of liquid through a Pasteur pipet.

Add 20λ R65 at 1/120 or R anti T ϕ 1/5; resuspend by above method

Incubate in ice for 20 min.

Resuspend after incubation by above method because cells settle during incubation.

Add 150λ D's PBS + 5% FCS; resuspend; centrifuge, flick, resuspend pellet.

Use 200λ D's PBS + 5% FCS to do 3 washes

Add 20λ R-goat anti rabbit at 1/20; incubate for 20 min and follow the above washing procedure.

Pellet, resuspend add 40λ FCS mix very well and smear for a thin smear
 20λ for a thick smear
or if stained at $1 \times 10^6/\text{well}$ add 100λ D's + 5% FCS to pellet, resuspend well and add to 1.9 ml D's + FCS ;
cells will be at $5 \times 10^5/\text{ml}$, but 0.2 ml in cytocentrifuge will have
 $\sim 10^5/\text{slide}$ cyto-centrifuge slide

Microtiter "V" plates - Cooke Engineering obtained through SP.

Microtiter Carriers - " "