

666-2117

MITHRAMYCIN STAINING PROTOCOL

used high power laser w. spectrum ---

1. Harvest drug-treated cells from spinner culture.
2. Centrifuge sample for 5 min at 1200 rpm.
3. Siphon off medium.
4. Resuspend cells in mithramycin solution using a bulb pipette and allow cells to stand 20 min at room temperature.
5. Analyze stained cells on FMF at 457 nm wavelength setting.

^{Wavelength, staining}
Mithramycin solution:

Mithramycin 100 $\mu\text{g/ml}$ in 25% ethanol and 15 mM MgCl_2 .

Stock @ 500 $\mu\text{g/ml}$
dilute 1 ml stock
1.25 ml ~~ETOH~~
2.75 ml PBS (without FCS)
Add 15 mg To 5 ml } $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$
3 mg/ml } MW 203.33