

HEPES BUFFER
for RPMI 1640
final conc. want 10m molar

F.W. 238.3

Make to 100X ∴

23.83 gm / 100ml in 25 gm bottle
105 ml

in d. H₂O & filter sterilize
pH to ~ 7.3

Use 1ml / 100ml RPMI

ORDER INFO:

HEPES BUFFER
GIBCO cat # 130440
25 gm bottle \$13.10
100 gm \$41.20

N-2-Hydroxyethylpiperazine-N'-
2-Ethanesulfonic acid